ACRON™ MC

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GC AMERICA INC.

ALSIPI, IL 60803 U.S.A.

GC AMERICA INC.
3737 West 127th Street, Alsip, Illinois 60803 U.S.A.

DISTRIBUTED by
GC EUROPE N.V.
Research-Park, Interleuvenlaan 13, B-3001 Leuven, Belgium

0500705

ACRON™ MC

MICROWAVE CURING DENTURE BASE RESIN

Prior to use, carefully read the instructions. For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS:

Construction of full and partial dentures.

CONTRAINDICATIONS:

Not to be used by or on individuals with known allergy to methacrylate. In case of allergy, refer to physician.

DIRECTIONS FOR USE:

- Coat the inside of the flask with thin film of Vaseline or similar lubricant.
- Place waxed denture into lower half of flask and fill it with 50/50 mix of plaster & stone.
- After the gypsum has hardened, apply Vaseline to the gypsum surface.
- Next fill up the top part of the flask with plaster/stone mix. Close the flask with the top deck.
 - NOTE:** To avoid plaster running into the bolt-holes, use the bolts to close the flask.
- After the gypsum has hardened, the wax is softened but be careful not to melt the wax. The softened wax should then be removed all in one piece. Next, the residual wax is thoroughly flushed out with boiling water containing an ordinary neutral detergent.
 - NOTE:** When utilizing the microwave flask, it is permissible to place in a microwave oven to initiate the wax removal process. Place on the turntable and irradiate for 1 minute on high. Then discard wax.
- After the wax has been removed, allow the gypsum surface to cool for about 5 minutes. This allows the steam to rise for better bonding of the separator. Next apply a separating medium (ex: GC ACRO-SEP or GC COE-SEP) to the gypsum surface and allow to dry.
- Mixing of powder and liquid for dough stage.
 - The standard powder/liquid ratio is 100g powder to 43mL liquid.
 - Place the required amount of liquid in the GC ACRON MC mixing jar.
 - Then place the corresponding amount of powder over the liquid while tapping the jar until excess liquid disappears.
 - NOTE:** Do not stir this mixture with a spatula to avoid air bubbles.
 - Close the lid of the mixing jar and leave for about 15-20 minutes at 74°F (23°C). At this point ACRON MC reaches a dough stage.
 - NOTE:** ACRON MC comes with a powder scoop and a 4mL liquid dropper graduated in 1mL increments. Average denture needs 2 scoops of powder and 8mL liquid.

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- Remove the dough in one piece from the mixing jar without producing air bubbles. Then, pack it into the model invested in the flask stored at room temperature. It is recommended not to touch the resin mixture directly because some individuals may be allergic to the liquid. Use a polyethylene sheet. Close the special microwavable flask. Trial closures are repeated 2 or 3 times until no flash is observed.
 - NOTE:** Trial packing is used to remove excess resin to avoid open vertical dimensions (= sur-occlusion).
- After trial packing, the flask should be fitted back together. Insert the bolts provided and apply pressure using the flask press. Tighten these bolts, under constant pressure, by using the supplied wrench.
- Cure ACRON MC with microwave irradiation for 3 minutes on high. For dentures with metal parts microwave for 1.5 minutes, then turn over for an additional 1.5 minutes. This is needed because microwaves cannot penetrate metal for curing the resin.
 - NOTE:** When using a 500 watt microwave oven use 100% power setting. When using a 1000 watt microwave oven use a 50% power setting. Failing to adjust the power setting may cause porosity in the denture.
- Cool the flask for more than 30 minutes at room temperature after irradiation, then place in cold water for complete cooling before deflasking (about 20 minutes).
- Follow conventional methods for deflasking the denture. This procedure is generally accomplished by the use of plaster-snips or a plaster knife.
- After deflasking put the denture into GC Roclean (not available in the U.S.A.), a plaster dissolvent, in order to totally clean the denture.
- Polish the denture by the usual method. There are several ways to obtain a high luster to the denture. GC does not recommend any specific brand of material for this purpose.

STORAGE:

Store in a cool dark place at room temperature (23°C/74°F) away from direct sunlight. (3 year shelf life guarantee).

PACKAGES:

Powder:	- refill	454 g	(1lb)	Liquid:	- refill	237 mL (8oz)
	- bulk	2.27 kg	(5lb)		- bulk	946 mL (32oz)

Microwavable flask.

Shades: Pink Veined, Pink Veined Light, Light Fibered +, Dark Pink Fibered, Clear

WARNING:

- Liquid is highly flammable. Avoid sources of ignition.
- Do not let the liquid contact your skin. When contacted, wash away thoroughly with soap and water. Some individuals may be allergic to the liquid.
- Irritation and sensitization have been reported in conjunction with products that contain methacrylates. To reduce risk of such occurrences, do not use this product with patients who have demonstrated sensitivity, avoid repeated and prolonged contact with uncured material, and discontinue use if sensitization occurs.

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