A-SILIKONBASERAT HYDROFILT AVTRYCKSMATERIAL

SÅ HÄR GÖR DU

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GC ASIA DENTAL PTE, LTD.

19 Loyang Way, #06-27 Singapore 508724
Tel: +65-6546-7588

HYDROPHILIC VINYL POLYSILOXANE IMPRESSION MATERIAL

HOW TO USE

Cavity Impression of Facing, Inlay, Onlay, etc.	Impression Technique	Tray Type	Recommended Consistencies
Impression of Crown or Bridge	Single Mix Single Impression	Stock or Custom Tray	Monophase
	Double Mix Single Impression		Monophase or Heavy + Injection or Regular
	Double Mix Double Impression		Putty + Injection or Regular
	Double Mix Single Impression		Monophase or Heavy + Injection or Regular
	Putty-Wash Impression		Putty + Injection or Regular
Impression of Partial Denture	Single Mix Single Impression	Stock or Custom Tray	Monophase
	Putty-Wash Impression	Stock Tray	Putty + Regular or Injection
Impression of Full Denture	Single Mix Single Impression	Stock or Custom Tray	Monophase


NEW CARTRIDGE LOADING AND DISPENSING

- Lift the release lever of the CARTRIDGE DISPENSER 2 (referred to as the dispenser hereafter) and pull the piston plunger all the way back into the dispenser. Lift the cartridge holder of the dispenser and load the cartridge, ensuring that the V-shaped notch on the flange of the cartridge is facing down. Push the cartridge holder down to hold the cartridge firmly in place.
- Remove the cartridge cap by rotating 1/4 turn anti-clockwise. Tilt the cap downward and peel it away from the cartridge. Gently squeeze the dispenser handle to extrude a small amount of material from the two openings at the end of cartridge. Make sure that base and catalyst come out evenly.
- Align the V-shaped notch on the rim of mixing tip to the V-shaped notch between the cartridge barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise to the end of the cartridge. The dispenser is now ready for use.
- Squeeze the handle several times to extrude the material. After use, do not remove the mixing tip as this will become the storage cap until next use. When replacing the mixing tip, rotate the collar on the mixing tip 1/4 turn anti-clockwise to align the V-shaped notch on the cartridge. Tilt the mixing tip downward and peel it away from the cartridge.
- Remove and replace the old mixing tip immediately prior to next use. Before attaching a new tip, gently extrude a small amount of material to ensure that base and catalyst are flowing evenly from both openings. If the materials should fail to extrude, remove any hardened materials from the end of the cartridge.
- To replace the cartridge, lift the release lever and retract the piston plunger fully. Remove the empty cartridge by lifting the cartridge holder and load a new one into the dispenser.

PHYSICAL PROPERTIES (Average)

Test Items	Injection Type 3: Light-Bodied Consistency	Regular Type 2: Med.-Bodied Consistency	Monophase Type 2: Med.-Bodied Consistency	Heavy Type 1: Heavy-Bodied Consistency
Total Working Time (Min., Sec.)	1'00"	1'15"	1'20"	1'15"
Setting Time (Min., Sec.)	1'45"	2'15"	2'15"	2'30"
Minimum Time in Mouth (Min., Sec.)	1'30"	1'30"	1'30"	1'30"
Recovery From Deformation (%)	99.7	99.5	99.3	99.3
Maximum Strain in Compression (%)	≤7.0	≤4.8	≤3.2	≤3.3
Linear Dimensional Change after 24 hrs (%)	≤0.2	≤0.2	≤0.2	≤0.2

Test conditions: Temperature (74°F/23°C ± 4°F/2°C) Relative humidity (50 ± 5%) (ISO 4823:1992(E))

The new-type cartridge must be used with the new GC CARTRIDGE DISPENSER 2.

DIRECTIONS FOR USE

- Attach an intraoral tip to the mixing tip or load material into a syringe.
- Inject the mixture onto the prepared teeth.
- Load the tray with the appropriate material (putty, monophase, heavy body or regular - depending on technique being used).
- Seat the loaded tray in the mouth within specified working times.
- Wait for set (1'30" minutes in the mouth).
- The obtained impression should be cleaned, then disinfected, utilizing a 2.5% or 3.4% glutaraldehyde, or other appropriate disinfectant, according to the manufacturers label recommendations.
- Remove the impression and pour a model immediately (if desired). Maximum time for pouring the model is 14 days.

NOTES

- When extruding EXAFAST NDS, care should be taken to avoid mixing or contact with the following materials. They may delay or prevent setting:
 - Catalyst for condensation type silicone impression materials
 - Polysulfide impression materials
 - Eugenol materials
 - Sulfur
 - Latex
 - Oil
 - Acrylates
 Also avoid moisture and glycerol when extruding.
- Fast set material may not be preferred for full arch multiple (more than 3) unit impressions where quick setting time could interfere with accuracy. Higher temperatures accelerate set; refrigeration retards set.
- Store in a dry and cool place.
- Avoid use of EXAFAST NDS with patients who have a history of hypersensitivity to silicone impression materials.
- EXAFAST NDS can be silver-or copper-plated.
- Care should be taken to avoid getting material on clothing. It is hard to remove.

HYDROPHILES VINYL-POLYSILOXAN ABFORMMATERIAL

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