Use only by a dental professional in the recommended indications. EN

GC Fuji PLUS™ EWT

RADIOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT WITH EXTENDED WORKING TIME

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Cementation of metal, ceramic and fiber posts

CONTRAINDICATIONS

- Direct pulp capping.
- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2,0/1,0
Mixing time (sec.)	20*
Working time (min., sec.) (at 23 °C) (from start of mixing)	3'30*
Final finishing contacting time (after sealing the restoration)	5'45*

- TOOTH PREPARATION**
 - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
 - Clean the prepared teeth with pumice and water. Note: Use GC Fuji PLUS CONDITIONER on the prepared tooth surface for 20 seconds for optimal adhesion.
 - Rinse thoroughly with water. Remove excess moisture by blotting with a cotton swab. Dry with an oil-free air. **DO NOT DESICCATE.** Prepared surface should appear moist (glistening).
- RESTORATION PREPARATION**
 - Prepare restorations as pretreated and handled according to the manufacturer’s instructions.

POWDER AND LIQUID DISPENSING

Use GC Fuji PLUS Liquid. Refer to the instructions for use of GC Fuji PLUS.

The standard powder to liquid ratio is 2.0g / 1.0g (1 level large scoop of powder to 3 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid).

- Do not dispense or powder, lightly tap the bottle against the cap. Do not shake or invert.
- Hold the liquid bottle vertically and squeeze gently.
- Dispense powder and liquid into the pad. Using the plastic spatula, add the liquid to the liquid. Mix rapidly for 20 seconds.
- MIXING**

Dispense powder and liquid onto the pad. Using the plastic spatula, add the liquid to the liquid. Mix rapidly for 20 seconds.

Note : When mixing larger amounts, divide the powder into 2 equal parts. Mix each portion with all of the liquid for 15 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-25 seconds (total 30-40 seconds).
- CEMENTATION**
 - Coat the internal surface of the restoration with sufficient cement and seal immediately. The working time is 3 minutes 30 seconds after cementation at 23 °C (73.4°F). Higher temperatures will shorten working time.
 - Maintain moderate excess cement when the excess cement feels rubbery.
 - Finishing can be started 5 minutes 45 seconds after sealing the restoration.

STORAGE

Recommended for optimal performance, store in a cool and dark place away from moisture and light.

SHADE
A3 (Vita® shade)
Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGE
Powder 15g with scoop

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or cosca butter can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT** mix powder or liquid with components of other glass ionomer cements.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS EWT in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS EWT from setting or bonding properly.
- This product is not indicated for filling or core build-up.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at:
http://www.georgelux.com
They can also be obtained from your supplier.

CLEANING AND DISINFECTING:
MULTI-STEP DELIVERY SYSTEMS. To avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.
DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects - Reporting:
You should be aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en
as well as to our internal vigilance system: vigilance@gc.dental
In this way you will contribute to improve the safety of this product.

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