GC Fuji PLUS®

RADIOOPAQUE REINFORCED GLASS IONOMER LUTING CEMENT

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Cementation of metal-based inlays, onlays, crowns and bridges
- Cementation of resin inlays, onlays, crowns and bridges
- Cementation of all ceramic inlays
- Cementation of high strength (zirconia based) all ceramic crowns and bridges
- Manufacturing of metal, ceramic and zircon posts

CONTRAINDICATIONS

- Pulp capping
- Avoid use of this product in patients with known allergies to glass ionomer cement, methacrylate monomer or methacrylate polymer.

DIRECTIONS FOR USE

Powder / Liquid ratio (g/g)	2.0/ 1.0
Mixing time (sec.)	20"
Working Time (23°C/73°F) (min., sec.)	2'30"
Finishing Time (37°C/99°F, after seating the restoration) (min., sec.)	4'30"

1) TOOTH PREPARATION

- Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
- Clean the prepared tooth with pumice and water.

Note:
Use GC Fuji PLUS CONDITIONER on the prepared tooth surface for 20 seconds for optimal adhesion.
- Rinse thoroughly with water. Remove excess moisture by blotting with a cotton pellet or gently blowing with an oil-free air. DO NOT DESICCATE. Prepared cavity should remain moist.
- RESTORATION PREPARATION**

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.
- POWDER AND LIQUID DISPENSING**

a) The standard powder to liquid ratio is 2.0 g / 1.0 g. (1 level large scoop of powder to 3 drops of liquid, or 1 level small scoop of powder to 1 drop of liquid)

b) For accurate dispensing of powder, lightly tap the bottle against the hand. Do not shake or invert.

c) Close bottles immediately after use.
- MIXING**

Dispense powder and liquid onto the pad. Using the plastic spatula, add all the powder to the liquid. Mix rapidly for 20 seconds.

Note:
When mixing larger amounts, divide the powder into 2 equal parts. Mix the first portion with all of the liquid for 15 seconds. Incorporate the remaining portion and mix the whole thoroughly for 15-25 seconds (total 30-40 seconds).
- CEMENTATION**

a) Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 30 seconds from start of mixing at 23°C (73°F). Higher temperature will shorten working time.

Note:
For more complex procedures requiring longer working time, it is recommended to mix powder and liquid on a cold mixing slab or use GC Fuji PLUS EWT.

b) Maintain moderate pressure.

c) Start removing excess cement when the excess cement feels rubbery.

d) Finishing can be started 4 minutes 30 seconds after seating the restoration.

STORAGE

Recommended for optimal performance, store in a cool and dark place (4-25°C (39.2-77.0°F)).

SHADES

Yellow, A3 (Vita® shade), Translucent Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

PACKAGES

- 1-1 package: 15 g powder, 8 g (7 mL) liquid, GC Fuji PLUS CONDITIONER 7 g (6.5 mL), a double-ended powder scoop, a plastic spatula, mixing pad (#23)
- Powder 15 g with scoop
- Liquid 8 g (7 mL)
- GC Fuji PLUS CONDITIONER 7 g (6.5 mL)

CAUTION

- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- DO NOT mix powder or liquid with components of other glass ionomer cements.
- GC Fuji PLUS CONDITIONER contains 2% Ferric Chloride. It is not recommended to use GC Fuji PLUS CONDITIONER in combination with products containing Epinephrine as this may stain oral tissue.
- Do not use GC Fuji PLUS in combination with eugenol containing materials as eugenol may hinder GC Fuji PLUS from setting or bonding properly.
- This product is not indicated for filling or core build-up.
 - Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.
- In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at: <http://www.gouropce.com> or <http://www.gamceria.com>

CLEANING AND DISINFECTING: MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged. DO NOT IMMERSE. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional / national guidelines.

Undesired effects- Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en, as well as to our internal vigilance system: vigilance@gc.dental. In this way you will contribute to improve the safety of this product.

Last revised: 11/2020

Prior to use, carefully read the instructions for use.

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