Prior to use, carefully read the instructions for use.

G-GOND™

ONE COMPONENT SELF-ETCHING LIGHT-CURED ADHESIVE

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- Bonding of light cured composites and acid modified composites (compomers) to denture dentures and metal frameworks.
- Bonding of dual cured luting and core build up composites to tooth structure as long as these materials are light cured.

CONTRAINDICATIONS

- Avoid use of this product in patients with known allergies to methacrylate monomer, methacrylate polymer, or acetone.

NOT TO BE USED

- In combination with chemically cured composite resin.
- In combination with eugenol containing materials as eugenol may hinder G-BOND from setting or bonding properly.
- In combination with desensitizers as desensitizers may hinder G-BOND from setting or bonding properly.
- In combination with dual cure luting and core build up composites in case these materials will not be light cured.

DIRECTIONS FOR USE**1. CAVITY PREPARATION**

Prepare tooth using standard techniques. Use rubber dam for protection.

Dry the prepared tooth surfaces by gently blowing with an air syringe.

Note: Do not use alcohol prep, use calcium hydroxide with.

2. APPLICATION

Direct to wetting, shake the bottle of G-BOND thoroughly (Fig. 1).

Dispense a few drops into an clean dispensing thru (Fig. 2).

Replace bottle cap immediately after the end of application using the DISPOSABLE APPLICATOR II (Fig. 4).

Apply G-BOND directly to the prepared enamel and dentin surfaces using the DISPOSABLE APPLICATOR II (Fig. 4).

Thick result should be a thin, rough, adhesive film with the appearance of frosted glass and which does not visibly move under further air pressure.

Light curing is not required.

When removing from refrigerator after long-term storage, leave the bottle of G-BOND to stand at room temperature for several minutes prior to use.

Rinse with water and gently dried before using a volatile solvent.

Apply G-BOND liberally as the material contains a volatile solvent (to be bonded) using a sponge or cotton pellet as residue material is hard to remove.

Should the applied material be contaminated with water, blood or saliva prior to light curing, wash and dry the tooth and repeat the procedure by using a fresh applicator.

3. LIGHT CURING

Light curing using a visible light curing unit (Fig. 7).

Irradiation time:

Halogen / LED (700 mW/cm²) 10 secondsHigh power LED (more than 1200 mW/cm²) 5 seconds

In case where the light guide tip is more than 10mm away from the surface to be irradiated, light cure for the following table:

Halogen / LED (700 mW/cm²) 20 secondsHigh power LED (more than 1200 mW/cm²) 10 seconds

Note: Light cure completely for an effective bonding strength.

Use a protective light shield or similar protective eye wear during light curing.

4. PLACEMENT OF LIGHT CURED COMPOSITES AND COMPOMERS

After having light cured the adhesive, follow manufacturer's instructions for placement of place, contour and light cure the composite resin material of choice.

4B. PLACEMENT OF DUAL CURED COMPOSITES

After having light cured the adhesive, make sure to light cure the dual cured composite material separately. Self curing only will result in improper adhesion.

5. FINISHING

Adjust restoration and polish using standard techniques.

STORAGE

Do not use for a prolonged period of time, store in refrigerator, otherwise store at room temperature (1-28°C) (33.8-82.4°F).

PACKAGE

1. Kit.

2. 5 mL liquid (1), dispensing dish (No.1) (1), DISPOSABLE APPLICATOR II (50)

3. 5 mL liquid (1)

Note: All packaging types may not be available in all countries.

CAUTION

G-BOND is flammable. Do not use near open flame. Keep away from sources of ignition. Do not store large quantities in one area. Keep away from direct sunlight.

G-BOND is volatile. Use in a well ventilated place.

In case of contact with eyes, flush immediately with water and seek medical attention.

In case of contact with oral tissue or skin, immediately remove with a sponge or cotton pellet. After the restorative treatment is finished, rinse thoroughly with water.

If the tissue contacted by the material turns white or forms a blister, advise the patient to leave the affected area undisturbed, until the mask disappears.

Avoid inhalation of dust. Avoid contact. It is recommended to apply COCOA BUTTER to the area where rubber dam cannot cover.

Avoid inhalation of digestion of material.

When spill on the table or floor, immediately wipe up with a dry cloth.

Do not mix with other products.

Keep away from children.

Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

In any case where the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a dentist.

Do not place under operating light.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <http://www.gcprouse.com>.

They can also be obtained from your supplier.

CLEANING AND DISINFECTING

If you become aware of any mid-level disinfectant, immediately after use inspect device and label for deterioration. Discard device if damaged.

DO NOT IMMERSER. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects - Reporting:

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:

https://ec.europa.eu/health/sensors/medical-devices/contacts_en

as well as to our customer care team, vigilance@gpc.dental

In this way you will contribute to improve the safety of this product.

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