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# G-CEM LinkAce™

SELF-ADHESIVE RESIN CEMENT

For use only by a dental professional in the recommended indications.

## RECOMMENDED INDICATIONS

- Cementation of all types of all ceramic, resin and metal-based inlays, onlays, crowns and bridges.
- Cementation of metal, ceramic, fiber posts, and cast post and cores.

## CONTRAINDICATIONS

- Pulp capping.
- Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

## DIRECTIONS FOR USE

A. CEMENTATION TECHNIQUE for ceramics, indirect composites and metal-based restorations

- TOOTH PREPARATION**
  - Prepare tooth in usual manner. For pulp capping, use calcium hydroxide.
  - Clean the prepared tooth with pumice and water. A clean surface is paramount for an optimal adhesion.
  - Rinse thoroughly with water. Dry thoroughly by gently blowing with oil free air. Prepared surface should be dry.
 

NOTE:

    - Do not use substances like H<sub>2</sub>O<sub>2</sub>, EDTA and/or disinfectants as these may lower the bond strength.
    - Special care should also be taken to remove imaging powders very thoroughly with pumice and water.
- RESTORATION PREPARATION**

Make sure that the restoration is pretreated and handled according to the manufacturer's instructions.
- DISPENSING**
  - Before dispensing the first time from a new syringe or after a long interval in between use, bleed the pastes from the syringe.
  - Remove the cap of the syringe by rotating 1/4 turn counterclockwise. Attach a GC Automix Tip by aligning the V shaped notch on the rim of the mixing tip with the V shaped notch between the syringe barrels. Push firmly to attach the mixing tip. Then rotate the colored collar of the mixing tip 1/4 turn clockwise.
  - Auxiliary material directly into the restoration.
 

NOTE:

The used mixing tip can be left on the syringe to serve as a storage cap until next use. For each application, place a new mixing tip.

- CEMENTATION**
    - Coat the internal surface of the restoration with sufficient cement and seat immediately. The working time is 2 minutes 45 seconds from start of mix at 23°C (73,4°F). Higher temperatures will shorten the working time.
 

NOTE:

Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation procedure.
  - Maintain moderate pressure until final polishing is done.
- EXCESS CEMENT REMOVAL**
  - Maintain moderate pressure and make sure the restoration remains in place.
  - To remove excess cement, light cure using a visible light curing unit for 1 second or allow the cement to self-cure until the excess cement feels rubbery.
  - Remove excess cement with appropriate instruments, holding the restoration in position.
 

NOTE:

Protect the margins of the restoration with a glycerin gel after removing the excess cement to prevent oxygen inhibition.
- FINAL SET**

While maintaining moderate pressure, light cure all surfaces / margins for 20 seconds each (Halogen/LED, intensity: 700 mW/cm<sup>2</sup>, covering wavelength: 430 to 480 nm). In case of the following restorations, let the material set for 4 minutes after removing excess cement.

  - Metal and metal coping restorations
  - Metal-free dark and opaque restorations
  - translucent metal-free restorations greater than 2mm in thickness
- FINAL POLISHING**

If needed, finish and polish margins with appropriate instruments.

- CEMENTATION TECHNIQUE for metal, ceramic, fiber posts, and cast post and cores**
  - TOOTH PREPARATION**

Prepare the post space according to manufacturer's instructions. The post space should be cleaned, rinsed and dried thoroughly using paper points. The post space should be dry.

NOTE:

    - Sodium hypochlorite solution is recommended to use for the chemical cleaning of the post space.
    - Do not use H<sub>2</sub>O<sub>2</sub> and/or EDTA to chemically clean the post space as these products may lower the bond strength.
  - POST PREPARATION**

Make sure that the post is pretreated and handled according to the manufacturer's instructions.
  - DISPENSING & CEMENTATION**
    - Before dispensing the first time from a new syringe or after a long interval in between use, bleed the pastes from the syringe.
    - Attach a GC Automix Tip for endo. For details, see A 3b) above.
    - Apply the extension tip to the mixing tip.
    - Extrude material directly into the post space.
    - Insert the post immediately into the post space within 1 minute after cement application.
 

NOTE:

      - The oral temperature accelerates the setting reaction of the cement.
      - Ambient light may shorten the working time. It is recommended to diminish the intensity or turn off your operation light during the cementation
      - Keep the tip immersed in the cement while dispensing to avoid trapping air.
      - DO NOT use lentulu spacers to place the cement in the post space as this may shorten the working time.
    - Maintain moderate pressure.
- EXCESS CEMENT REMOVAL**

Continue to maintain moderate pressure making sure the post remains in place, and remove excess cement with appropriate instruments. In case of a cast post and core, tack cure each all surface for 1 second each with a visible light curing unit.
- FINAL SET**
  - While maintaining moderate pressure, light cure all surfaces/margins for 20 seconds each (Halogen/LED, intensity: 700 mW/cm<sup>2</sup>, covering wavelength: 430 to 480 nm).
  - Let the material set for another 4 minutes.

**STORAGE**  
Recommended for optimal performance, store at room temperature away from heat, moisture and direct sunlight (4-25°C) (39,2-77,0°F).

**SHADE**  
A2 (Vita® shade), A03 (opaque), B01 (opaque) and Translucend Vita® is a registered trademark of Vita Zahnfabrik, Bad Säckingen, Germany.

- PACKAGE**
- G-CEM LinkAce syringe 4.6 g (2.7 mL) (2), GC Automix Tip Regular (15), GC Automix Tip for endo with extension tip (5)
  - GC Automix Tip Regular (20)
  - GC Automix Tip Regular (10)
  - GC Automix Tip for endo with extension tip (10)
- CAUTION**
- In case of contact with oral tissue or skin, remove immediately with a sponge or cotton soaked in alcohol. Flush with water. To avoid contact, a rubber dam and/or COCOA BUTTER can be used to isolate the operation field from the skin or oral tissue.
  - In case of contact with eyes, flush immediately with water and seek medical attention.
  - Do not use G-CEM LinkAce in combination with eugenol containing materials as eugenol may hinder G-CEM LinkAce from setting or bonding properly.
  - Use a protective light shield or similar protective eye wear during light curing.
  - This product is not indicated for filling or core building.
  - In rare cases the product may cause sensitivity in some people. If any such reactions are experienced, discontinue the use of the product and refer to a physician.
  - Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheets available at <http://www.gceurope.com> or for the Americas <http://www.gcamerica.com>

They can also be obtained from your supplier.

**CLEANING AND DISINFECTING**  
MULTI-USE DELIVERY SYSTEMS: to avoid cross-contamination between patients this device requires mid-level disinfection. Immediately after use inspect device and label for deterioration. Discard device if damaged.

**DO NOT IMMERS**E. Thoroughly clean device to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Undesired effects- Reporting:  
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: [https://ec.europa.eu/growth/sectors/medical-devices/contacts\\_en](https://ec.europa.eu/growth/sectors/medical-devices/contacts_en) as well as to our internal vigilance system: [vigilance@gcc.dental](mailto:vigilance@gcc.dental) in this way you will contribute to improve the safety of this product.

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