GC Acrylic Primer
Light-curing Primer for bonding composites to acrylics

For use only by a dental professional in the recommended indications.

INTENDED USE

Increases the adhesiveness of light-curing composites to conventional acrylic resins used in dental laboratory procedures (e.g. modification of denture teeth / denture base resins).

CONTRAINDICATIONS

Avoid use of this product in patients with known allergies to methyl methacrylate monomer, poly methyl methacrylate polymer or acetone.

DIRECTIONS FOR USE

- Preparation of adhesion area:
Roughen the surface by sandblasting with aluminum oxide (50µm at 2bar) or by roughening by tungsten burs. Clean the surface by steam cleaner or clean with oil free air. The prepared surface should not be touched anymore.
- Applying GC Acrylic Primer:
Dispose of a few drops of GC Acrylic Primer into a dispensing dish.
Wet the adhesion area with GC Acrylic Primer by using a brush and light-cure the area to polymerize the material. Irradiation time and curing unit

Curing Unit	LABOLIGHT DUO	STEPLIGHT SL-I	LABOLIGHT LV-III, LV-II
	Step-mode	Full-mode	Pre-cure
GC A. Pr.	20sec.	3min.	20sec.

NOTE:

- Light-cure completely for an effective bonding strength.
- Use a protective light shield or similar protective eye wear during light-curing.
- Use fresh material and avoid leaving the material dispensed for a long period of time. Its high volatilization may cause insufficient curing or may decrease the bonding strength if it has been left too long on a dispensing dish.
- Remove excess material to avoid any cause of decrease of the bonding strength.

- Layering of composite resin materials.
Before layering make sure that the to be layered surface is free of any impurities. Place light-curing composite resin materials on the adhesion area and polymerize them in accordance with the manufacturer's instruction for use (e.g. GRADIA PLUS (GC)).

STORAGE
Recommended for optimal performance, store in a cool and dark place (-4°-25°C / 39.2-77.0 °F) away from high temperatures or direct sunlight.

PACKAGES
GC Acrylic Primer, Bottle 6mL

CAUTION

- In rare cases, the product may cause sensitivity in some people. If such reaction occurs, discontinue the use of the product and refer to a physician.
- The device cannot be used by other than by a dental professional.
- The device cannot be used for any indication other than in accordance with the manufacturer's instruction for use.
- When cutting and polishing the prosthesis, wear a safety goggles, use a dust collector and wear a dust mask to avoid inhalation of cutting dust.
- Don't reuse a brush with other products.
- When carrying out the work, to provide sufficient ventilation to prevent headache etc. If it were inhaled by human, move to fresh location.
- Avoid inhalation/ ingestion of material.
- Do not use near naked flame. Keep away from sources of ignition. Keep away from direct sunlight.
- If not in use for a prolonged period of time, store in refrigerator, otherwise store at room temperature. GC Acrylic Primer is to be used at room temperature.
- GC Acrylic Primer is volatile. Use in a well-ventilated place. Apply to adhesion area and close the bottle quickly after application.
- Before closing the liquid bottle we advise to clean the nozzle. Always vertically close the bottle.
- In case of contact with eyes, rinse immediately with water and seek medical attention.
- In case of contact with oral tissue or skin, immediately remove with a sponge or cotton soaked in alcohol. Rinse immediately with water.
- In case of surface spills, immediately wipe off with a dry cloth.
- Do not mix with other similar products.
- Dispose all waste in accordance with local regulations.
- Personal protective equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the safety data sheet available at: <http://www.gceurope.com>. They can also be obtained from your supplier.

Last revised: 06/2018

CE 0086

MANUFACTURED BY
GC DENTAL PRODUCTS CORP.
2-285 Toriguchi-cho, Kasugai,
Aichi 486-0844, Japan

DISTRIBUTED BY
GC CORPORATION
7-1 Hasunuma-cho, Itabashi-ku,
Tokyo 174-8585, Japan

GC EUROPE N.V.
Researchpark Haasrode-Leuven 1240,
Interleuvenlaan 33,
B-3001 Leuven, Belgium TEL: +32 16 74 10 00

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- Layering of composite resin materials. Before layering make sure that the to be layered surface is free of any impurities. Place light-curing composite resin materials on the adhesion area and polymerize them in accordance with the manufacturer's instruction for use (e.g. GRADIA PLUS (GC)).

STORAGE

Recommended for optimal performance, store in a cool and dark place (-4°-25°C / 39.2-77.0 °F) away from high temperatures or direct sunlight.

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