# Master all your challenges

*initial*<sup>TM</sup>

MC Classic Line

- EN Technical Manual
- DE Technisches Handbuch
- FR Manuel Technique
- IT Manuale tecnico
- ES Manual Técnico
- PL Manual Técnico

JULY 2021

**GC**

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
EN Thank you for choosing our GC Initial MC veneering ceramic.  
GC Initial MC is a veneering ceramic with EN ISO 6872 for the ceramic veneering of precious and non-precious metal substructures with a coefficient of thermal expansion of the metal substructure in the range of  $13,8 - 14,9 \cdot 10^{-6} \text{ K}^{-1}$  (25-500°C).  
Also suitable for electroforming copings and for the sintering technique to create inlays, onlays and veneers on refractory dies.  
This instruction for use will give you an impression of how easy it is to get a convincing aesthetic result with a minimum workload that highlights the excellent features of this veneering ceramic.

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---

## EN Safe and reliable

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- Reliable 3-layering technique
- Fast, easy and reliable matching of all Vita-shades
- Secure and stable firing behaviour
- Minimal shrinkage
- Outstanding handling properties
- Compatible with all GC Initial MC Expert & IQ powders
- Colour cross-match with all GC Initial veneering ceramics

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## EN Physical facts

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- GC Initial MC is a veneering ceramic with a leucite reinforced microstructure
- Highly homogeneous
- High flexural strength
- Gingiva friendly
- Antagonist friendly
- Easy mechanical polishing in the mouth
- Firing stable opalescence
- Natural light behaviour

Dental ceramic complying with EN ISO 6872 for the ceramic veneering of precious and non-precious metal substructures with a coefficient of thermal expansion of the metal substructure in the range of  $13,8 - 14,9 \cdot 10^{-6} \text{ K}^{-1}$  (25-500°C).

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---

## EN Safety Notes

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- GC Initial MC is a veneering ceramic with a lucite reinforced microstructure
- Inhaling ceramic dust can cause damage to the lungs.
- The ceramic powder must be stored in a tightly sealed container and used according to the manufacturer's instructions.
- Always wear a mask with a P1 filter or use an adequate and fully suction unit.

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## EN Before you start ...

## DE Bevor Sie starten ...

- GC Initial MC is a veneering ceramic with a lucite reinforced microstructure
- Make sure that the CTE of the alloy you use stays strictly within the range of  $13,8 - 14,9 \cdot 10^{-6} \text{ K}^{-1}$  (25-500°C).
- Calibrate your furnace according to the manufacturer's instruction in order to obtain the best results and to be able to make optimal use of the properties.
- As each dental ceramic is sensitive to contamination, take care that you work in a clean environment.

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## DE Komponenten des GC Initial MC Classic Sets

### EN INmetalbond

Intermediate coating between the alloy and the ceramic to block escaping metal oxides and as CTE Buffer.

#### **Paste Opaque**

To individualise and to cover the framework according to the desired Vita shade.

#### **Dentin**

To build up the core and to define the basic colour.

#### **Enamel**

Used in the upper third of the restoration to match the enamel.

#### **Glaze Powder / Liquid**

To give the restoration a natural glaze.

#### **Modelling Liquid**

To mix the ceramic materials during build-up.

#### **Opacus Dentin**

Opacus Dentin is used to intensify the opacity and to increase the chroma.

#### **CLW (Clear Window)**

To add transparent effects.

#### **Correction Powder**

For small corrections after the glaze firing.

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EN When preparing the metal substructure please follow the alloy manufacturer's instructions for use. In case, try to look up the latest information at the manufacturer's homepage or contact your local sales representative!

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- EN Important note when using non-precious alloys:
- Please note that non-precious alloys are in general more sensitive to their manipulation. Please check the compatibility before use. Best use the alloys recommended by GC.
  - Pd free dental alloys („bio“ type) containing zinc (Zn) should be sandblasted, oxidized. After firing soak in a clean, warm acid bath for approx. 5- 6 minutes.
  - To prepare the metal surface, please use sharp cross-cut drills.
  - Grind the surface with slow speed and low pressure.
  - Grind the surface in one direction only to avoid overlaying layers (to avoid air bubbles after the first Opaque firing).

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## EN INmetalbond

## DE INmetalbond

EN The bonder is used as an intermediate coat between the alloy and the ceramic.  
It blocks escaping metal oxides and neutralizes differences in the expansion coefficient (CTE).  
INmetalbond can be used on all alloys within the given CTE range.

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## EN Opaque bake

## DE Opakerbrand

EN Whether you first applied INmetalbond (recommended) or directly the Paste Opaque onto the metal structure, you always should follow the following steps:  
Apply a thin layer of Wash Opaque (semi-masking) evenly to the coping or bridge framework. Please use a Paste Opaque brush.  
Ensure that this first layer is not too thick. If a too thick layer of Paste Opaque is applied, uncontrollable amounts of liquid evaporate during the heating stage and can cause the occurrence of bubbles.  
After the Wash Opaque bake is finalized, the second Opaque bake is carried out.  
Do not dry the Paste Opaque too quickly. If the recommended preheating and firing times are not followed, small cavities can open up at the surface or the Opaque may lift off from the substructure. This happens, if the Opaque liquid evaporates too quickly by an improper or accelerated firing schedule!

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## EN Dentin & Incisal bake

### DE Dentin- & Schneidebrand

EN Mix the Dentin and Incisal powders with the Modelling Liquid. Add Modelling Liquid until the ceramic meets your desired consistency.

Start with the Dentin core and continue the build-up with Dentin and Enamel. In order to compensate the firing shrinkage the build-up needs to be slightly oversized. In the same way, corrections can be done with your preferred Dentin and Enamel shades.

Note: Never mix GC Initial MC products with ceramics from other suppliers or other GC Initial products that are not explicitly pointed out to be compatible with GC Initial MC.

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## EN Stain & Glaze bake

## DE Malfarben- & Glasurbrand

EN For glazing and individualisations, only use the recommended GC Initial materials.

GC recommends:

1. GC Initial Spectrum Stains and Glaze Liquid,  
or
2. GC Initial IQ, Lustre Pastes ONE / NF Gum

Note 1: Stains and Glaze Liquid can only be mixed with the corresponding special liquids. Apply the Stains with a pointed brush.

Note 2: The GC Initial IQ, Lustre Pastes ONE / NF Gum can be used as well for characterisation. Please refer to the IFU included in the set or download from the GC Europe web page: <https://europe.gc.dental/nl-BE/products/initialiqlustrepastesone>

Note

The degree of gloss is determined by the firing temperature and the holding time. Increasing the firing temperature / holding time will lead to a glossier surface! Lowering the temperature / holding time will create a matte appearance.

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## EN Optional Steps

## DE Optionale Schritte

EN In case you need to intensify the chroma of the desired shade, you can add an Opaqus Dentin in order to give the restoration an increased intensity.

To compensate the light reflection at the margin of the sub-structure you should apply the Opaqus Dentin at the indicated area (image 2).

You can use an Opaqus Dentin as the first layer when you start the build-up process.

Please take into account that excessive use of Opaqus Dentin will lead to a high chroma and an increased opacity. To get acquainted with the ratio of Opaqus Dentin and Dentin we recommend to make a trial beforehand.

### CLW (Clear Window)

In case certain areas require translucency, CLW (Clear Window) can be applied or mixed with the other ceramic. Also here GC recommends to make a trial beforehand to get acquainted with the effect of CL-W.

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# EN Firing Instructions

## DE Brenntabellen

### INITIAL MC GENERAL FIRING INSTRUCTION

	START TEMP. STARTTEMPERATUR	DRYING TIME VORTROCKNUNGSZEIT	HEAT RATE AUFHEIZRATE	RELEASE VACUUM RELEASE VACUUM	FIRING TEMP. BRENNTEMPERATUR	HOLDING TIME HALTEZEIT	VACUUM [HPA] VAKUUM [HPA]	EXTENDED COOLING LANGZEITABKÜHLUNG
INmetalbond	550°C	6 min	80°C/min	Yes Ja	980°C	1 min	50	-
Wash Opaque* Wash-Brand*	550°C	6 min	80°C/min	Yes Ja	940°C	1 min	50	-
Paste Opaque Pastenopaker	550°C	6 min	80°C/min	Yes Ja	930°C	1 min	50	-
1 <sup>st</sup> Dentin Firing 1. Dentinbrand	580°C	6 min	55°C/min	Yes Ja	890°C	1 min	50	-
2 <sup>nd</sup> Dentin Firing 2. Dentinbrand	580°C	6 min	55°C/min	Yes Ja	880°C	1 min	50	-
Glaze Firing (without Stain/Glaze) Glanzbrand (ohne Malfarben / Glasurmasse)	600°C	2 min	55°C/min	No Nein	890°C	1 min	-	-
Glaze Firing (with Stain/Glaze) Glanzbrand (mit Malfarben / Glasurmasse)	480°C	2 min	45°C/min	No Nein	850°C	1 min	-	-
Correction Powder Firing Korrekturmasse	450°C	4 min	45°C/min	Yes Ja	770°C	1 min	50	-

#### EN Important note:

\* In case non-precious alloys are used without INmetalbond, the firing temperatures of the Wash Opaque firing should be increased with 20°C. This is not applicable in case INmetalbond is applied.

#### Note:

The given firing temperatures are recommendations only and may differ depending on the furnace. The temperatures have to be adjusted to the furnace you are using and according to the result after firing. Always make sure your furnace is correctly calibrated.

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EN Troubleshooting  
DE Fehlerbehebung

FAULT	POSSIBLE CAUSE	CORRECTIVE ACTION
Ceramic opaque is too porous	<ul style="list-style-type: none"> <li>• preheating temperature is too high</li> <li>• heat up is too fast</li> <li>• vacuum pump starts too late</li> <li>• vacuum level reached is too low</li> </ul>	<ul style="list-style-type: none"> <li>• reduce preheating temperature</li> <li>• reduce heat up speed</li> <li>• correct vacuum starting time</li> <li>• check vacuum pump and furnace for leak tightness and function</li> </ul>
Color too light, opaque	<ul style="list-style-type: none"> <li>• preheating temperature is too high</li> <li>• vacuum level too low</li> <li>• firing temperature too low</li> </ul>	<ul style="list-style-type: none"> <li>• reduce preheating temperature</li> <li>• correct vacuum level</li> <li>• check and adapt firing temperature</li> </ul>
Ceramic surface is too smooth / too rough	<ul style="list-style-type: none"> <li>• end temperature is too high /too low and (or) holding time too long /too short</li> </ul>	<ul style="list-style-type: none"> <li>• decrease/increase the end temperature and (or) decrease/increase holding time</li> </ul>
Ceramic surface not shiny	<ul style="list-style-type: none"> <li>• end temperature is too low and (or) holding time is too short</li> </ul>	<ul style="list-style-type: none"> <li>• increase end temperature and (or) increase holding time</li> </ul>
Edges and contours too round	<ul style="list-style-type: none"> <li>• end temperature too high and (or) holding time too long</li> </ul>	<ul style="list-style-type: none"> <li>• reduce end temperature and (or) decrease holding time</li> </ul>
Cracks / pressure stress cracks: horizontal cracks in the incisal area or pontics	<ul style="list-style-type: none"> <li>• wrong alloy (CTE) used</li> <li>• wrong cooling down</li> <li>• wrong shape of metal substructure</li> </ul>	<ul style="list-style-type: none"> <li>• check CTE of alloy</li> <li>• follow cooling down recommendations according to CTE</li> <li>• modeling reduced anatomical tooth shape as metal substructure</li> </ul>
Cracks / tension stress cracks: vertical crack (craquelure cracks)	<ul style="list-style-type: none"> <li>• sharp edges in the metal substructure</li> <li>• wrong alloy (CTE)</li> <li>• wrong cooling down</li> </ul>	<ul style="list-style-type: none"> <li>• avoid sharp edges in the metal substructure</li> <li>• check the CTE of the alloy</li> <li>• adapt cooling down according to CTE, see recommendations</li> </ul>

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GC NORDIC AB - Finnish Branch
GC NORDIC AB - GC Nordic Danish Branch
GC NORDIC AB
GC AUSTRIA GmbH - Swiss Office
GC AUSTRALASIA DENTAL PTY LTD
GC ASIA DENTAL PTE. LTD.