EN Prior to use, carefully read the instructions for use.


GC KALORE™ LIGHT-CURED RADIOPAQUE UNIVERSAL COMPOSITE RESTORATIVE

For use only by a dental professional in the recommended indications.


- RECOMMENDED INDICATIONS
1. Direct restorative for Class I, II, III, IV, V cavities.
2. Direct restorative for wedge-shaped defects and root surface cavities.
3. Direct restorative for veneers and diastema closure.

- CONTRAINDICATIONS
1. Pulp capping.
2. Avoid use of this product in patients with known allergies to methacrylate monomer or methacrylate polymer.

- DIRECTIONS FOR USE
1. Shade Selection
Clean the tooth with pumice and water. Shade selection should be made prior to isolation.
2. Application
For bonding GC KALORE to enamel and / or dentin, use a light-cured bonding system such as G-BOND, G-animal Bond or G-Premio BOND (Fig. 1).

- Clinical Hints
a. In the case of small cavities, Restoring using a one shade technique. In most cases the use of one Universal shade alone will be sufficient.
b. In the case of large and/or deep cavities, In most cases a multi shade layering technique will give the best aesthetic results.

GC KALORE shade combination chart for multiple layers in deep and/or large cavities. Table with columns for shade (A1-A4, B1-B3, C1-C3, D1-D2, WT, DT, NT) and rows for Opaque, Translucent, and Universal.

- For details of shades, refer to the following section of SHADES.
5. Contouring before Light Curing
6. Light Curing
Light cure GC KALORE using a light curing unit (Fig. 3). Keep light guide as close as possible to the tooth surface.

GC KALORE Irradiation Time and Effective Depth of Cure. Table with columns for irradiation time (10, 20, 30 sec) and effective depth (3.0, 3.5, 2.5, 2.0, 1.5 mm).

- 1. Material should be placed and light cured in layers.
2. Lower light intensity may cause insufficient curing or discoloration of the material.
7. Finishing and Polishing
Finishing and polish using diamond burs, polishing points and discs.
SHADES
25 shades
15 Universal shades (color code on unit/cap / syringe label: green)

- PACKAGES
I. Unittips
1. Refill
a. Pack of 20 tips (each in 11 shades) (0.16mL per tip)
b. Pack of 10 tips (each in 15 shades) (0.16mL per tip)

- CAUTION
1. In case of contact with oral tissue or skin, remove immediately with cotton or a sponge soaked in alcohol.
2. In case of contact with eyes, flush immediately with water and seek medical attention.

- CLEANING AND DISINFECTING
MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients this device requires a high level of disinfection.
2. Do not use GC KALORE in combination with eugenol containing materials as eugenol may hinder GC KALORE from setting.

- UNDESIRABLE EFFECTS - REPORTING
If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link:
https://ec.europa.eu/growth/sectors/medical-devices/contacts_en

GC DENTAL PRODUCTS CORP. 2-285 Tomimatsu-cho, Kasugai, Aichi 486-0844, Japan. RESPONSIBLE MANUFACTURER IN CANADA: GC AMERICA INC. 3737 West 127th Street, Aliso Viejo, CA 92683 U.S.A. MANUFACTURED BY GC CORPORATION 7-1-1 Hasegawa-cho, Itabashi-ku, Tokyo 174-8585, Japan. EU: GC EUROPE N.V. Researchpark Haardtsele-Leuven 1240, Interleuvenlaan 33, B-3001 Leuven, Belgium. GC AMERICA INC. 3737 West 127th Street, Aliso Viejo, CA 92683 U.S.A. GC SOUTH AMERICA Rua Heliópolis, 399, Santana - São Paulo, SP, BRASIL. CE 2797

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