PATTERN RESIN™ LS

ACRYLIC RESIN FOR PATTERNS

For use only by a dental professional in the recommended indications.

RECOMMENDED INDICATIONS

- To make the pattern of:
 - tapered coping for telescopic crowns and other telescope crowns.
 - inlay, palatal and lingual bar, connector and clasp.
 - Maryland bridge.
- Temporary splint for soldering.
- In the mouth for core build up pattern.
- Construction of implant superstructures.

CONTRAINDICATIONS

In susceptible individuals, sensitisation to the product cannot be excluded. Use of the product should be discontinued if allergic reactions are observed.

WARNING

Contains methyl methacrylate. When used directly in the mouth for core build up particular care should be taken to avoid contact with the skin or oral mucosa. In case of accidental contact, immediately remove the offending material and flush the surrounding area liberally with water.

DIRECTIONS FOR USE

1. Powder and liquid preparation

Dispense adequate amounts of powder and liquid into each rubber cup (Fig. 1).

2. Building the coping (for secondary crown patterns of tapered crowns or telescopic crowns)

- Moisten the tip of the brush with the liquid and pick up a small amount of powder (Fig. 2).
- Deposit this mixture on the crown and progressively build up the secondary crown pattern by repeating the procedure (Fig. 3).
 - In this case, no separating medium is required, but when working on a stone model, GC Acro Sep, GC MULTI-SEP or Vaseline should be used as separator.
 - Clean the brush with liquid after use or during the work if necessary.

3. Checking the internal surface of the pattern

After the resin has set (3 min. at 23°C/73°F), remove the pattern carefully from the internal crown and check the margin and the internal surface (Fig. 4).

4. Fine adjustment

- Return the pattern to the internal crown and make a uniform thickness (0.3-0.4 mm) by using a bur (Fig. 5).
- Remove the pattern and check the thickness with the measuring device (Fig. 6).

5. Waxing veneer and cut out

- Perform the wax-up with inlay wax (Fig. 7).

Note: There is no need to clean the surface of the pattern resin with a cleaning agent before waxing.
- Cut out the veneer for the facing (Fig. 8).
- Place the retention beads (Fig. 9).

6. Investing, Burning out and Casting

For a single crown, all procedures can be performed in the usual manner. For a larger pattern, hold at 250°C (482°F) for 1 hour before raising to the final burn out temperature (Fig. 10).

STORAGE

Recommended for optimal performance, store in a cool place away from direct sunlight.

NOTES

- Do not mix with other chemically cured resins.
- Do not return the residual powder and liquid to the bottle.
- In case of contact with skin wash off with soap and water.

PACKAGES

- Powder: 100g, 1kg
- Liquid: 100g (105mL), 250g (262mL)
- 1-1 Pkg: 100g powder, 100g (105mL) liquid with accessories

CAUTION

- The liquid is highly FLAMMABLE. Avoid sources of ignition.
- In case of contact with eyes, flush immediately with water and seek medical attention.
- Personal Protective Equipment (PPE) such as gloves, face masks and safety eyewear should always be worn.

CLEANING AND DISINFECTING RECOMMENDATION

MULTI-USE DELIVERY SYSTEMS: To avoid cross-contamination between patients the bottles and measuring devices require mid-level disinfection. Immediately after use inspect the bottles, measuring devices and label for deterioration. Discard if damaged. DO NOT IMMERSSE. Thoroughly clean bottles and measuring devices to prevent drying and accumulation of contaminants. Disinfect with a mid-level registered healthcare-grade infection control product according to regional/national guidelines.

Some products referenced in the present IFU may be classified as hazardous according to GHS. Always familiarize yourself with the Safety Data Sheets available at:

<http://www.gceurope.com> or for the Americas: <http://www.gcamerica.com>

They can also be obtained from your supplier.

Undesired effects - Reporting

If you become aware of any kind of undesired effect, reaction or similar events experienced by use of this product, including those not listed in this instruction for use, please report them directly through the relevant vigilance system, by selecting the proper authority of your country accessible through the following link: https://ec.europa.eu/growth/sectors/medical-devices/contacts_en as well as to our internal vigilance system: vigilance@gc.dental

In this way you will contribute to improve the safety of this product.

Last revised: 08/2019

GC 2797 **Rx Only**

MANUFACTURED BY
GC AMERICA INC.
3127 West 127th Street, Atlap, IL 60803 U.S.A.
TEL: +1-708-587-0900
www.gcamerica.com

RESPONSIBLE MANUFACTURER IN CANADA
GC AMERICA INC.
3127 West 127th Street, Atlap, IL 60803 U.S.A.

DISTRIBUTED BY
GC CORPORATION
76-1 Hasunuma-cho, Itabashi-Ku, Tokyo 174-8585, Japan

GC EUROPE N.V.
Biosciencepark Hazardode-Louven 1340, Interleuvenlaan 33,
B-3001 Louven, Belgium TEL: +32-16-74-10-00

GC SOUTH AMERICA
Rua Heliodora, 300, Santana - São Paulo, Sp - Brasil
CEP: 02022-051 - TEL: +55-11-2625-9965
CNPJ: 08.279.999/0001-61
RESP. TEC. Mayara De Santes Ribeiro - CRO/SP: 105.982

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