# GC VEST-G

JULY 2019

NR3605  
950130  
30000072

- EN** Phosphate-bonded investment  
for all dental alloys
- DE** Phosphatgebundene Einbettmasse  
für Dentallegierungen
- FR** Revêtement au phosphate pour  
tous les alliages dentaires
- IT** Rivestimento fosfato per qualsiasi  
lega dentale
- ES** Revestimiento al fosfato para  
todo tipo de aleaciones dentales

**GC**

# GC VEST-G

EN

## - PHOSPHATE-BONDED INVESTMENT FOR ALL DENTAL ALLOYS -

GC VEST-G is an improved phosphate-bonded investment for all dental alloys, for crowns and bridges. The predosed pack and controllable expansion make it easy to have more precise casting which modern technique requires.

### ADVANTAGES


- The expansion is adjustable to shrinkage of different casting alloys.
- Fine powder provides smooth and dense surface of the cast.
- Easy to deflask.
- Enough strength is kept and no harmful effect appears on casting alloys at high temperature.

### PHYSICAL PROPERTIES (at 23°C, Liquid concentration: 100%)

L/P ratio	24 ml/100g
Setting time	14 minutes
Total expansion	3.26%
Compressive strength (dry)	12MPa

### LIQUID CONCENTRATION

Kind of alloy	liquid : distilled water	Total expansion
Non-precious metal	85-100% : 15-0%	3.16-3.26%
Semi-precious metal	50% : 50%	2.64%
Precious metal	33% : 67%	2.40%

### DIRECTIONS FOR USE

#### (1) Preparation

Line a ring with a dry non-asbestos type liner (1mm thickness) such as CASTING RING LINER.


If necessary, dilute the GC VEST-G liquid with distilled water depending on the kind of alloy and work.

We highly recommend to use a wetting agent for wax patterns.

#### (2) Mixing

Powder	Liquid
60g	14.4ml
150g	36.0ml


#### Liquid measuring

14.4 ml for 60g  
36.0 ml for 150g

Mix powder and liquid for 60 seconds with a vacuum-mixer (420 rpm).

#### Heating

(depending on the ringsize and quantity of rings)

(3) *Drying and wax elimination*

1 hour keeping at room temperature		
1 hour heating up slowly to 250°C (→ 4°C/Min)		
1 hour keeping at 250°C		
1 hour heating up slowly to final temperature (→ 8°C/Min)		
1 hour keeping at final temperature		
Final temperature:	700-750°C for	precious metal
	800-850°C for	precious metal for porcelain, semi-precious metal, non-precious metal.

**NOTES**

- Clean the bowl, spatula,... very well after use.
- Before use, keep the powder and liquid at room temperature (about 23°C).
- Liquid may freeze when exposed to temperature below 0°C. Once frozen, the liquid cannot be used anymore.
- Ring method is recommended.

**SAFETY RECOMMENDATIONS AND HAZARD WARNINGS**

1. Investment materials contain quartz.  
Do not inhale dust!  
Risk of lung damage (Silicosis, lung cancer).  
Advise: wear protection mask Type FFP 2 – EN 149:2001.  
Open the investment material bag with scissors and avoid the formation of dust when filling into the mixing bowl.  
Rinse the empty investment material bag with water before disposal.
2. Remove dust from your working place only when it is wet.
3. To avoid the formation of dust when removing the investment material from the casting ring, place the cooled casting ring into water for a short time.
4. When sandblasting the cast object, always use a fine-dust filter extraction system.

**PACKAGES**

Powder: 6kg Carton (60g x 100)  
6kg Carton (150g x 40)  
12kg Tin

Liquid: 900 ml Bottle

Expiry date: 3 years from the manufacturing date.

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**TAISEI DENTAL MFG. Co, LTD.**