

One shade fits all



Essentia[®]
Universal shade
from GC

Universal resin composite



GC

Essentia[®] Universal shade from GC


One shade fits all

With **Essentia Universal shade**, GC throws the door to simplification wide open. Essentia Universal exhibits an excellent chameleon effect and **seamlessly blends with the surrounding tooth structure, regardless of the shade of the cavity.**


Available in **three viscosities** - all strong enough to be used **in all cavity classes** - Essentia Universal is the ultimate solution for posterior restorations.

Thanks to a proprietary filler treatment and to the use of ultra-fine particles, **the physical properties of Essentia HiFlo & Essentia LoFlo are even higher** than those of Essentia Universal (paste type), confirming that all products can be safely used **for all indications.**

Essentia Universal Shade (paste) blending in teeth of different shades...



Flexural strength & 3-body wear of the three viscosities of Essentia Universal shade (Source: GCC R&D data)



Preparation of the small Class I cavities



Restoration using Essentia HiFlo Universal



Post-operative view showing great blending with natural teeth



Air abrasion of the non-carious lesions



Restoration with **Essentia LoFlo Universal**




Post-operative view after polishing



Essentia HiFlo is a fluid composite with a **high wettability**, perfect for easy adaptation to preparations and flow into narrow cavities. It offers great handling properties together with superior strength and a high radiopacity of over 300% facilitating the diagnosis of secondary caries.



Essentia LoFlo is an injectable composite with a **low flowability**, optimal when a thixotropic material is needed – for instance for cervical restorations. It promises minimal wear, ensuring long-lasting restorations even on the occlusal surface.



Essentia is a paste composite with a simplified 7-shade concept. Its Universal shade displays a great blending effect together with a **packable consistency**: perfect for mono-shade posterior restorations.



A direct approach was selected for this challenging post-endodontic case



Restoration with **everX Posterior** as a dentine core and **Essentia Universal** as the outer shell



Six-month recall showing good shade integration into the surrounding tissues

Essentia® Universal shade from GC

One shade fits all



Essentia HiFlo

Universal light-cured composite
with high flowability



Essentia LoFlo

Universal light-cured composite
with low flowability



Essentia

Universal light-cured composite
with a packable consistency

900739 Essentia HiFlo Syringe, U (1 x 2mL)
900741 Essentia LoFlo Syringe, U (1 x 2mL)
900970 Essentia Syringe, U (1 x 2mL)
900998 Essentia Unitips, U (15 x 0.16mL)

Also discover G-Premio BOND

GC's universal bonding featuring zero compromises

012687 G-Premio BOND Kit 5mL
012694 G-Premio BOND Bottle Refill 5mL
012688 G-Premio BOND 3-Bottle Pack
009034 G-Premio BOND Unit Doses, 50 pieces (0.1mL)



e s s e n t i a

GC EUROPE N.V.

Head Office
Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info.gce@gc.dental
<http://www.gceurope.com>

GC UNITED KINGDOM Ltd.

Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info.uk@gc.dental
<http://uk.gceurope.com>

GC