

<b>VITA® Classic</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A4</b>	<b>B1</b>	<b>B2</b>	<b>B3</b>	<b>B4</b>	<b>C1</b>	<b>C2</b>	<b>C3</b>	<b>C4</b>	<b>D2</b>	<b>D3</b>	<b>D4</b>
Opak	AO1	AO2	AO3	AO3	AO3	AO1	AO2	AO3	AO3	AO2	AO3	AO3	AO3	AO1	AO2	AO3
<b>CORE</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A4</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>	<b>A3.5</b>	<b>A2</b>	<b>A3</b>	<b>A3</b>	<b>A3.5</b>	<b>A1</b>	<b>A2</b>	<b>A3</b>
Enamel	JE	JE	AE	AE	AE	JE	JE	AE	AE	JE	AE	AE	AE	JE	JE	AE

### Multi-Farb-Technik-Korrespondenz

G-ænial® **A'CHORD**

Enamel



Opak



CORE



Enamel



CORE



Enamel



**Spritze**  
(1 x 2.1ml/4.0g) 10006739  
**Unitips**  
(15 x 0.16ml/0.33g) 10006753



**CORE-Farben** A1  
**VITA® Classic** A1 B1 D2

10006740

10006754



A2  
A2 B2 C1 D3

10006741

10006755



A3  
A3 B3 C2 C3 D4

10006742

10006756



A3.5  
A3.5 B4 C4

10006743

10006757



A4  
A4

## G-ænial® A'CHORD

### Opak-Farben

BOW



10006750

AO1



10006747

AO2



10006748

AO3



10006749

### Standard-Farben

BW



10006746

A5



10006744

A6



10006745

### Enamel-Farben

JE



10006751

AE



10006752

**Spritze**  
(1 x 2.1ml/4.0g)  
**Unitips**  
(15 x 0.16ml/0.33g)

10006764

10006761

10006762

10006763

10006760

10006758

10006759

10006765

10006766