

VITA® Classic	A1	A2	A3	A3.5	A4	B1	B2	B3	B4	C1	C2	C3	C4	D2	D3	D4
Opaco	AO1	AO2	AO3	AO3	AO3	AO1	AO2	AO3	AO3	AO2	AO3	AO3	AO3	AO1	AO2	AO3
Base (CORE)	A1	A2	A3	A3.5	A4	A1	A2	A3	A3.5	A2	A3	A3	A3.5	A1	A2	A3
Esmalte	JE	JE	AE	AE	AE	JE	JE	AE	AE	JE	AE	AE	AE	JE	JE	AE

Correspondencia de la técnica estratificada

G-ænial® **A'CHORD**

Esmalte

Opaco

CORE

Esmalte

CORE

Esmalte



Jeringa
(1 x 2.1ml/4.0g)
Unitips
(15 x 0.16ml/0.33g)

10006739

10006753



Colores BASE
VITA® Classic
A1 B1 D2

10006740

10006754



A2
A2 B2 C1 D3

10006741

10006755



A3
A3 B3 C2 C3 D4

10006742

10006756



A3.5
A3.5 B4 C4

10006743

10006757



A4
A4

G-ænial® A'CHORD

Colores opacos

BOW



10006750

AO1



10006747

AO2



10006748

AO3



10006749

Colores estándar

BW



10006746

A5



10006744

A6



10006745

Colores de esmalte

JE



10006751

AE



10006752

Jeringa
(1 x 2.1ml/4.0g)
Unitips
(15 x 0.16ml/0.33g)

10006764

10006761

10006762

10006763

10006760

10006758

10006759

10006765

10006766