

G-Premio BOND

from GC



Make it **simple**

Same bottle, same procedure
for all your adhesive work



Direct
composite



Luting



Core
build-up



Repair

GC

G-Premio BOND from GC


One-component light-cured **universal adhesive** for

- ✓ Direct bonding
- ✓ Repair cases
- ✓ Indirect luting
- ✓ Core build-up
- ✓ Hypersensitivity treatment


that can be used in **all etching modes**




Zero compromises on performance - **Zero** compromises on handling



3 functional monomers:
high bond strength to dentine,
enamel & various substrates



Low film thickness:
decreased risk of discolouration
& perfect marginal fit



Excellent penetration:
low post-operative sensitivity



Handy flip-top Cap!



Reduced waste
thin nozzle -
up to 300 drops
per bottle



G-Premio BOND DCA
for dual-curing
with GRADIA CORE
& post luting

WEP Code	EEP Code	Description
012687	012689	G-Premio BOND 5mL Kit
012694	012695	G-Premio BOND 5mL Bottle Refill
012688	012690	G-Premio BOND 3-bottle Pack
009034	009035	G-Premio BOND Unit Dose 50 pcs
009552		G-Premio BOND DCA 3mL Bottle



GC EUROPE N.V.
Head Office
Researchpark
Haasrode-Leuven 1240
Interleuvenlaan 33
B-3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info.gce@gc.dental
http://www.gceurope.com

GC UNITED KINGDOM Ltd.
Coopers Court
Newport Pagnell
UK-Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info.uk@gc.dental
http://uk.gceurope.com

