



Gradia[®] Direct LoFlo from GC.

High Viscosity Flowable
Microfilled Hybrid
Composite with
HDR technology.



Gradia[®] Direct LoFlo, the **first flowable** to be considered as a **real restorative**, not only for class V or as a liner, but also for class I, class II and class III restorations.

Introducing the new HDR, High-Density Radiopaque, technology, GC presents a new flowable restorative with unsurpassed physical properties.

stability – reliability – beauty


GC




High wear resistance and low shrinkage for your safety

Gradia® Direct LoFlo achieves excellent stability thanks to the unique High-Density Radiopaque (HDR) Technology. **Wear resistance and fracture toughness** can easily keep up with modern posterior composites and, on top of that, offer up to **55% less shrinkage** than the leading flowable composite on the market.

% Volume Polymerization Shrinkage



Wear (µm)



Reliability

The HDR Pre-polymerized Filler with Nano-silica filler technology significantly increase **durability, polishability and radiopacity**. The incorporation of fluoro-aluminosilicate glass brings you protection against secondary decay in addition.



Beautiful handling

Gradia® Direct LoFlo **wettability** allows for perfect **adaptability**, and adapts easily without slumping while moving easily with any type of instrument or brush of your choice. Due to its unique composition, similar to Gradia® Direct, the 7 available shades of Gradia® Direct LoFlo blend invisibly into the surrounding tooth structure.



Picture from Prof. Myasaki

Packaging

2 x 1,3 g syringes, 4 metal needle tips and 2 light protective caps.

Available in 7 shades: A1, A2, A3, A3.5, A03 (Opaque Dentin), CV (Cervical Shade), and BW (Bleach White).

Accessories

Dispensing tips, metal needle type or plastic type, 20 pcs.



GC EUROPE N.V.

Head Office
Interleuvenlaan 33
B - 3001 Leuven
Tel. +32.16.74.10.00
Fax. +32.16.40.48.32
info@gceurope.com
www.gceurope.com

GC UNITED KINGDOM Ltd.

12-15, Coopers Court
Newport Pagnell
UK - Bucks. MK16 8JS
Tel. +44.1908.218.999
Fax. +44.1908.218.900
info@uk.gceurope.com
www.uk.gceurope.com

